


DIATHERMY/ BIPOLAR

Specialized electrosurgical instruments for precise cutting, dissection, and coagulation. Includes monopolar and bipolar forceps, electrodes, handpieces, and connecting cables. Bipolar instruments deliver focused hemostasis with minimal thermal spread for delicate tissue; monopolar options provide versatility for broader dissection. High-grade insulation, ergonomic grips, and reliable conductivity ensure safe, consistent performance and exceptional surgeon control.



Bipolar Forceps

Straight



Fine tip

LS-DIA/034	101mm long
LS-DIA/035	127mm long
LS-DIA/036	152mm long
LS-DIA/037	178mm long
LS-DIA/038	185mm long


Medium tip

LS-DIA/039	101mm long
LS-DIA/040	127mm long
LS-DIA/041	152mm long
LS-DIA/042	178mm long
LS-DIA/043	185mm long

Heavy tip

LS-DIA/044	101mm long
LS-DIA/045	127mm long
LS-DIA/046	152mm long
LS-DIA/047	178mm long
LS-DIA/048	185mm long


Adson Bipolar Forceps



LS-DIA/003a	152mm long
-------------	------------

Bipolar Forceps

Angled



Fine tip

LS-DIA/049	101mm long
LS-DIA/050	127mm long
LS-DIA/051	152mm long
LS-DIA/052	178mm long
LS-DIA/053	185mm long

Medium tip

LS-DIA/054	101mm long
LS-DIA/055	127mm long
LS-DIA/056	152mm long
LS-DIA/057	178mm long
LS-DIA/058	185mm long

Heavy tip

LS-DIA/059	101mm long
LS-DIA/060	127mm long
LS-DIA/061	152mm long
LS-DIA/062	178mm long
LS-DIA/063	185mm long

Monel Handle




LS-DIA/028	145mm long
------------	------------

Needle Handles



LS-DIA/027

Abbey Diathermy Needle



LS-DIA/078


Finger Switch Pencil

With spatula and 3 metre lead




LS-DIA/068

Monel Blade



LS-DIA/029	Fully insulated
LS-DIA/030	Non-insulated

Needles



Straight

LS-DIA/064	Stout
LS-DIA/064a	Medium

Angled

LS-DIA/064b	Stout
LS-DIA/064c	Medium

Diathermy Arm

With ball



LS-DIA/069	Curved
LS-DIA/069a	Straight

Spatula Blade

For diathermy finger

LS-DIA/067a



Needle Point


For diathermy finger

LS-DIA/067b



Turner Warwick Diathermy Scissors Riches Artery Forceps

Screw joint, curved



Straight


LS-DIA/072	190mm long
LS-DIA/073	230mm long
LS-DIA/074	250mm long

Curved

LS-DIA/075	190mm long
LS-DIA/076	230mm long
LS-DIA/077	250mm long

LS-DIA/014	127mm long
LS-DIA/015	152mm long
LS-DIA/016	178mm long
LS-DIA/017	185mm long
LS-DIA/018	228mm long

Wilson Hey Artery Forceps




Straight

LS-DIA/023	152mm long
LS-DIA/024	178mm long
LS-DIA/025	185mm long
LS-DIA/026	228mm long

Curved


LS-DIA/019	152mm long
LS-DIA/020	178mm long
LS-DIA/021	185mm long
LS-DIA/022	228mm long

Adson Forceps



LS-DIA/003	152mm long
------------	------------

Dissecting Forceps



Serrated

LS-DIA/004	152mm long
LS-DIA/005	178mm long
LS-DIA/006	203mm long
LS-DIA/007	228mm long
LS-DIA/008	254mm long

1-2 teeth

LS-DIA/009	152mm long
LS-DIA/010	178mm long
LS-DIA/011	203mm long
LS-DIA/012	228mm long
LS-DIA/013	254mm long