


# GENERAL INSTRUMENTS

The essential foundation of every surgical toolkit, our General Instruments category covers the versatile, multi-purpose tools used across virtually all specialties and procedures. This collection includes fundamental instruments for tissue handling, hemostasis, cutting, dissection, suturing, and exposure—such as forceps, scissors, needle holders, retractors, and clamps. Designed for reliability, precision, and everyday durability, these instruments support the core tasks of surgery from initial incision to final closure.




Adson Artery Forceps

Box joint, 180mm long



Serrated jaws

|            |          |
|------------|----------|
| LS-GI/034  | Straight |
| LS-GI/034a | Curved   |




1-2 teeth

|            |          |
|------------|----------|
| LS-GI/034b | Straight |
| LS-GI/034c | Curved   |

Cairns Artery Forceps

Box joint, half serrated jaws  
145mm long



Serrated jaws

|           |          |
|-----------|----------|
| LS-GI/415 | Straight |
|-----------|----------|




1-2 teeth

|           |                |
|-----------|----------------|
| LS-GI/020 | Curved         |
| LS-GI/021 | Curved to side |

Crile Artery Forceps

Box joint, serrated jaws



Straight

|            |            |
|------------|------------|
| LS-GI/014  | 140mm long |
| LS-GI/014a | 160mm long |




Curved


|            |            |
|------------|------------|
| LS-GI/015  | 140mm long |
| LS-GI/015a | 160mm long |

Cushing Artery Forceps

Box joint, fine half serrated jaws, 140mm long




|           |          |
|-----------|----------|
| LS-GI/001 | Straight |
| LS-GI/002 | Curved   |



Dunhill Artery Forceps

Box joint, fine serrations



Straight


|            |            |
|------------|------------|
| LS-GI/065  | 135mm long |
| LS-GI/065a | 190mm long |

Curved

|           |            |
|-----------|------------|
| GI/043    | 135mm long |
| LS-GI/468 | 180mm long |

Lawson Tait Artery Forcep


Box joint, straight jaw



|              |            |
|--------------|------------|
| LS-GI/SH5924 | 140mm long |
|--------------|------------|

Halstead Mosquito Forceps

Standard pattern, haemostatic, box joint



Serrated jaw

|           |                         |
|-----------|-------------------------|
| LS-GI/007 | Straight 125mm long     |
| LS-GI/008 | Curved 125mm long       |
| LS-GI/009 | Fully curved 125mm long |

1-2 teeth


|           |                     |
|-----------|---------------------|
| LS-GI/010 | Straight 125mm long |
| LS-GI/012 | Straight 190mm long |
| LS-GI/011 | Curved 125mm long   |
| LS-GI/013 | Curved 190mm long   |

Micro, 100mm long

|           |                     |
|-----------|---------------------|
| -GI/003   | Straight, serrated  |
| LS-GI/004 | Curved, serrated    |
| LS-GI/005 | Straight, 1-2 teeth |
| LS-GI/006 | Curved, 1-2 teeth   |


Heiss Artery Forceps

Box joint, 205mm long



Serrated jaws

|           |          |
|-----------|----------|
| -GI/016   | Straight |
| LS-GI/017 | Curved   |



1-2 teeth

|          |          |
|----------|----------|
| -GI/016a | Straight |
| GI/017b  | Curved   |

Kelly Artery Forceps

Box joint, half serrated jaws



Straight


|           |            |
|-----------|------------|
| GI/022    | 140mm long |
| LS-GI/023 | 160mm long |



Curved

|           |            |
|-----------|------------|
| -GI/024   | 140mm long |
| LS-GI/025 | 160mm long |

Fraser Kelly Artery Forceps Curved




Straight

|            |       |
|------------|-------|
| LS-GI/468a | 180mm |
|------------|-------|


Kocher Artery Forceps

Box joint, 1-2 teeth



Straight

|           |            |
|-----------|------------|
| LS-GI/052 | 130mm long |
| LS-GI/053 | 150mm long |
| LS-GI/054 | 180mm long |
| LS-GI/055 | 200mm long |
| LS-GI/056 | 230mm long |




Curved

|           |            |
|-----------|------------|
| LS-GI/057 | 150mm long |
| LS-GI/058 | 180mm long |
| LS-GI/059 | 200mm long |
| LS-GI/060 | 230mm long |

Moynihan Artery Forceps

Box joint, half-serrated jaws



Straight


|           |            |
|-----------|------------|
| LS-GI/018 | 145mm long |
| LS-GI/068 | 180mm long |

Curved

|           |            |
|-----------|------------|
| LS-GI/019 | 145mm long |
| LS-GI/069 | 180mm long |

Rochester-Pean Artery Forceps

Box joint




Straight

|           |            |
|-----------|------------|
| -GI/044   | 165mm long |
| LS-GI/045 | 190mm long |
| LS-GI/046 | 205mm long |
| LS-GI/047 | 240mm long |

Curved

|           |            |
|-----------|------------|
| LS-GI/048 | 165mm long |
| LS-GI/049 | 180mm long |
| LS-GI/050 | 205mm long |
| LS-GI/051 | 240mm long |

Dieffenbach Bulldog Clamp



Straight, serrated

|            |                           |
|------------|---------------------------|
| LS-GI/112  | 1.2mm fine jaw, 40mm long |
| LS-GI/113  | 1.2mm fine jaw, 45mm long |
| LS-GI/108  | 1.8mm jaw, 30mm long      |
| LS-GI/109  | 1.8mm jaw, 35mm long      |
| LS-GI/110  | 2.5mm jaw, 40mm long      |
| LS-GI/111  | 2.5mm jaw, 45mm long      |
| LS-GI/107  | 6mm heavy jaw, 60mm long  |
| LS-GI/455  | 20mm jaw, 75mm long       |
| LS-GI/455a | 30mm jaw, 85mm long       |
| LS-GI/455b | 45mm jaw, 105mm long      |
| LS-GI/455c | 65mm jaw, 120mm long      |

Straight, plain

|           |                           |
|-----------|---------------------------|
| -GI/129   | 1.2mm fine jaw, 40mm long |
| LS-GI/131 | 1.2mm fine jaw, 45mm long |
| LS-GI/121 | 1.8mm jaw, 30mm long      |
| LS-GI/123 | 1.8mm jaw, 35mm long      |
| LS-GI/125 | 2.5mm jaw, 40mm long      |
| LS-GI/127 | 2.5mm jaw, 45mm long      |

Curved, serrated



|            |                           |
|------------|---------------------------|
| -GI/119    | 1.2mm fine jaw, 40mm long |
| LS-GI/120  | 1.2mm fine jaw, 45mm long |
| LS-GI/115  | 1.8mm jaw, 30mm long      |
| LS-GI/116  | 1.8mm jaw, 35mm long      |
| LS-GI/117  | 2.5mm jaw, 40mm long      |
| LS-GI/118  | 2.5mm jaw, 45mm long      |
| LS-GI/114  | 6mm heavy jaw, 60mm long  |
| LS-GI/456  | 20mm jaw, 75mm long       |
| LS-GI/456a | 30mm jaw, 85mm long       |
| LS-GI/456b | 45mm jaw, 105mm long      |
| LS-GI/456c | 65mm jaw, 120mm long      |

Curved, plain

|           |                           |
|-----------|---------------------------|
| -GI/130   | 1.2mm fine jaw, 40mm long |
| LS-GI/132 | 1.2mm fine jaw, 45mm long |
| LS-GI/122 | 1.8mm jaw, 30mm long      |
| LS-GI/124 | 1.8mm jaw, 35mm long      |
| LS-GI/126 | 2.5mm jaw, 40mm long      |
| LS-GI/128 | 2.5mm jaw, 45mm long      |

Roberts Artery Forceps


Screwjoint, curved on flat, serrated



|               |            |
|---------------|------------|
| S C- G I /070 | 230mm long |
|---------------|------------|

Spencer Wells Artery Forceps

Box joint, serrated jaws



Straight


|           |            |
|-----------|------------|
| -GI/035   | 127mm long |
| LS-GI/036 | 152mm long |
| LS-GI/037 | 177mm long |
| LS-GI/038 | 203mm long |

Curved

|           |            |
|-----------|------------|
| LSGI/039  | 127mm long |
| LS-GI/040 | 152mm long |
| LS-GI/041 | 177mm long |
| LS-GI/042 | 203mm long |

OTHER SIZES AVAILABLE

Blalock Bulldog Clamp



Serrated, with hole in handle

Straight

|           |                     |
|-----------|---------------------|
| LSGI/480  | 11mm jaw, 38mm long |
| LS-GI/482 | 14mm jaw, 48mm long |
| LS-GI/493 | 17mm jaw, 60mm long |
| LS-GI/481 | 27mm jaw, 70mm long |

Curved

|           |                     |
|-----------|---------------------|
| LSGI/478  | 11mm jaw, 38mm long |
| LS-GI/479 | 14mm jaw, 48mm long |
| LS-GI/491 | 17mm jaw, 60mm long |
| LS-GI/492 | 27mm jaw, 70mm long |

Plain, with hole in handle

Straight


|           |                     |
|-----------|---------------------|
| LS-GI/486 | 11mm jaw, 38mm long |
| LS-GI/487 | 14mm jaw, 48mm long |
| LS-GI/488 | 17mm jaw, 60mm long |
| LS-GI/489 | 27mm jaw, 70mm long |

Curved

|           |                     |
|-----------|---------------------|
| LS-GI/483 | 11mm jaw, 38mm long |
| LS-GI/484 | 14mm jaw, 48mm long |
| LS-GI/485 | 17mm jaw, 60mm long |
| LS-GI/490 | 27mm jaw, 70mm long |

Dandy Artery Forceps


Box joint, curved in plane



|               |            |
|---------------|------------|
| S C- G I /076 | 140mm long |
|---------------|------------|

Carmalt Angiotribe Forceps

Box joint, atraumatic, light pattern




Straight

|           |            |
|-----------|------------|
| LSGI/062  | 165mm long |
| LS-GI/063 | 200mm long |
| LS-GI/064 | 230mm long |

Curved

|           |       |
|-----------|-------|
| LS-GI/061 | 165MM |
|-----------|-------|


Adsons Dissecting Forceps



127mm long

|           |                              |
|-----------|------------------------------|
| LS-GI/257 | 1-2 teeth                    |
| LS-GI/258 | Lightweight, micro           |
| LS-GI/259 | Serrated                     |
| LS-GI/260 | Serrated, lightweight, micro |


Block End Dissecting Forceps



Serrated

|            |            |
|------------|------------|
| LSGI/255   | 127mm long |
| LS-GI/256  | 152mm long |
| LS-GI/256a | 203mm long |

Bonney Dissecting Forceps




178mm long

|           |                     |
|-----------|---------------------|
| LSGI/267  | Serrated, 1-2 teeth |
| LS-GI/268 | 1-2 teeth           |

Charnley-McIndoe Dissecting Forceps


150mm long



|             |           |
|-------------|-----------|
| LS-OR/N045a | 2-3 teeth |
|-------------|-----------|

DeBakey Atraumatic Dissecting Forceps

Straight, atraumatic jaws



1.5mm

|            |            |
|------------|------------|
| LSGI/238a  | 150mm long |
| LS-GI/240a | 240mm long |

2.5mm

|            |            |
|------------|------------|
| LS-GI/238b | 150mm long |
| LS-GI/240b | 240mm long |


2.8mm

|           |            |
|-----------|------------|
| GI/238    | 150mm long |
| LS-GI/241 | 180mm long |

|            |            |
|------------|------------|
| LS-GI/238c | 150mm long |
| LS-GI/240c | 240mm long |
| LS-GI/241c | 180mm long |

Dissecting Forceps

Turn over And, serrated jaws




|           |            |
|-----------|------------|
| LS-GI/242 | 127mm long |
| LS-GI/243 | 152mm long |

|           |            |
|-----------|------------|
| LS-GI/244 | 178mm long |
| LS-GI/245 | 203mm long |

Dunhill Dissecting Forceps


Light, fine serrated jaws



|           |            |
|-----------|------------|
| LS-GI/457 | 145mm long |
|-----------|------------|

Fredrick Coole Dissecting Forceps

Extra fine serrated jaws




125mm long

|           |          |
|-----------|----------|
| LSGI/066  | Straight |
| LS-GI/067 | Curved   |

Gillies Dissecting Forceps


1-2 teeth



|           |            |
|-----------|------------|
| LS-PL/052 | 152mm long |
|-----------|------------|

Jackson Burrows Dissecting Forceps


140mmlong



|             |                              |
|-------------|------------------------------|
| LS-OR/N037  | Tongue and groove, 1-2 teeth |
| LS-OR/N037a | Tongue and groove only       |

Jean Dissecting Forceps


Serrated jaws,140mmlong



|           |                       |
|-----------|-----------------------|
| LSGI/269  | Fine                  |
| LS-GI/270 | 1-2 teeth             |
| LS-GI/271 | 2-3 teeth, 178mm long |

Lanes Dissecting Forceps

1-2 teeth on serrated jaws, 140mm long




|           |            |
|-----------|------------|
| LS-GI/246 | 127mm long |
| LS-GI/247 | 152mm long |
| LS-GI/248 | 178mm long |

|           |                       |
|-----------|-----------------------|
| LSGI/249  | 2-3 teeth, 178mm long |
| LS-GI/250 | 2-3 teeth 203mm long  |

McIndoe Dissecting Forceps


152mm long



|           |           |
|-----------|-----------|
| LSGI/476  | Serrated  |
| LS-GI/477 | 1-2 teeth |


McIndoe Modified Dissecting Forceps

152mmlong, 1mm tip,modifiedforfinepowerfuldissection with strong action tip



|         |          |
|---------|----------|
| SH/5961 | Serrated |
|---------|----------|

Treves Dissecting Forceps




1-2 teeth

|           |            |
|-----------|------------|
| LS-GI/251 | 127mm long |
| LS-GI/252 | 152mm long |

2-3 teeth

|           |            |
|-----------|------------|
| LSGI/253  | 127mm long |
| LS-GI/254 | 152mm long |

Waugh Dissecting Forceps



Fine, serrated


|           |            |
|-----------|------------|
| LS-GI/261 | 152mm long |
| LS-GI/262 | 178mm long |
| LS-GI/263 | 203mm long |

Fine, 1-2 teeth

|           |            |
|-----------|------------|
| LSGI/264  | 152mm long |
| LS-GI/265 | 178mm long |
| LS-GI/266 | 203mm long |

Allis Tissue Forceps

Box joint



150mm long


|           |           |
|-----------|-----------|
| LSGI/272  | 3-4 teeth |
| LS-GI/273 | 4-5 teeth |

190mm long

|           |           |
|-----------|-----------|
| LS-GI/414 | 5-6 teeth |
|-----------|-----------|

Littlewood Tissue Forceps


Box joint, 2-3 teeth



|           |            |
|-----------|------------|
| LS-GI/275 | 190mm long |
|-----------|------------|

Babcock Tissue Forceps


Box joint, serrated jaws



|           |            |
|-----------|------------|
| LSGI/276  | 160mm long |
| LS-GI/277 | 200mm long |
| LS-GI/278 | 230mm long |

Duval Tissue Forceps


Box joint



|           |                       |
|-----------|-----------------------|
| LS-GI/279 | 13mm wide, 203mm long |
|-----------|-----------------------|

Lane Tissue Forceps


Box joint, 1-2 teeth, fenestrated



|           |            |
|-----------|------------|
| LS-GI/280 | 127mm long |
| LS-GI/281 | 152mm long |
| LS-GI/282 | 197mm long |

Rutherford Morrison Tissue Forceps


Box joint



|           |                       |
|-----------|-----------------------|
| LS-GI/283 | 4-5 teeth, 152mm long |
|-----------|-----------------------|

Shaw Tissue Forceps


Box joint



|           |                       |
|-----------|-----------------------|
| LS-GI/458 | 1-2 teeth, 178mm long |
|-----------|-----------------------|

Guys Gallabin Suture-Holding Forceps


Screw joint, angled flat blades



|            |            |
|------------|------------|
| LS-GI/450  | 125mm long |
| LS-GI/4402 | 150mm long |

Russian Dissecting Forceps


With fenestrated end



|                |            |
|----------------|------------|
| LS-GI/SH-5503b | 150mm long |
| LS-GI/SH-5503  | 200mm long |
| LS-GI/403      | 230mm long |
| LS-GI/SH-5503a | 250mm long |

Rampley Sponge Holding Forceps


Box joint, straight, serrated jaws



|           |            |
|-----------|------------|
| LSGI/353  | 178mm long |
| LS-GI/354 | 241mm long |

Foerster Sponge Holding Forceps

Box joint, straight, serrated jaws



|            |            |
|------------|------------|
| LS-GI/354a | 178mm long |
|------------|------------|


Harrison Bowl Steriliser Forceps

Screw joint, angled, double jaw

|           |            |
|-----------|------------|
| LSGI/365  | 356mm long |
| LS-GI/366 | 457mm long |

Lister Sinus Forceps


Box joint, serrated tip




|           |            |
|-----------|------------|
| LSGI/367  | 127mm long |
| LS-GI/368 | 152mm long |
| LS-GI/369 | 178mm long |

Cheate Sterilizer Forceps

Screw joint




|           |            |
|-----------|------------|
| LS-GI/364 | 267mm long |
|-----------|------------|



Hunter Splinter Forceps

Serrated grooved pointed jaws



|           |            |
|-----------|------------|
| LS-GI/360 | 115mm long |
|-----------|------------|

Michel Skin Holding Forceps


With magazine for Michel clips



|           |                       |
|-----------|-----------------------|
| LS-GI/358 | 1-2 teeth, 127mm long |
|-----------|-----------------------|

Childe Approximating Forceps


With magazine for Michel clips



|           |            |
|-----------|------------|
| LS-GI/469 | 180mm long |
|-----------|------------|

Michel Forceps


Box joint, double-ended, for inserting and removing clips stainless steel



|           |            |
|-----------|------------|
| LS-GI/357 | 130mm long |
|-----------|------------|

Michel Clip Extracting Forceps


Box joint



|           |            |
|-----------|------------|
| LS-GI/359 | 125mm long |
|-----------|------------|


Robins Anchoring Forceps

140 mm long



|            |                                     |
|------------|-------------------------------------|
| LS-GI/453  | 5mm with holder for diathermy leads |
| LS-GI/453a | 13mm with holder for suction tubing |

Dressing Scissors



Blunt/Blunt  
Straight

|           |            |
|-----------|------------|
| LS-GI/133 | 130mm long |
| LS-GI/136 | 146mm long |
| LS-GI/139 | 165mm long |

Curved

|           |            |
|-----------|------------|
| LS-GI/142 | 130mm long |
| LS-GI/143 | 146mm long |
| LS-GI/144 | 165mm long |
| LS-GI/145 | 185mm long |
| LS-GI/404 | 147mm long |

With one serrated blade

Blunt/Sharp  
Straight

|           |            |
|-----------|------------|
| LS-GI/134 | 130mm long |
| LS-GI/137 | 146mm long |
| LS-GI/140 | 165mm long |

Sharp/Sharp  
Straight

|           |            |
|-----------|------------|
| LSGI/135  | 130mm long |
| LS-GI/138 | 146mm long |
| LS-GI/141 | 165mm long |

Curved


|           |            |
|-----------|------------|
| LS-GI/146 | 130mm long |
|-----------|------------|

Dressing Scissors, Left-Handed

With one black bow,145mm long

Blunt/Blunt

|           |          |
|-----------|----------|
| LSGI/192  | Straight |
| LS-GI/193 | Curved   |



Blunt/Sharp

|           |          |
|-----------|----------|
| LSGI/194  | Straight |
| LS-GI/195 | Curved   |

Sharp/Sharp

|           |          |
|-----------|----------|
| LSGI/196  | Straight |
| LS-GI/197 | Curved   |

Littler Dissecting Scissors


Fine pattern,with eye for suture



|            |            |
|------------|------------|
| LS-PL/001f | 115mm long |
|------------|------------|

Kilner Dissecting Scissors

Fine pattern



|            |            |
|------------|------------|
| LS-PL/001a | 115mm long |
|------------|------------|

Metzenbaum Scissors



Straight

|            |            |
|------------|------------|
| LSGI/405   | 180mm long |
| LS-GI/406  | 147mm long |
| LS-GI/406a | 200mm long |
| LS-GI/406b | 230mm long |
| LS-GI/406c | 250mm long |
| LS-GI/406d | 280mm long |

Curved

|            |            |
|------------|------------|
| LS-GI/407  | 185mm long |
| LS-GI/408  | 147mm long |
| LS-GI/408a | 200mm long |
| LS-GI/408b | 230mm long |
| LS-GI/408c | 250mm long |
| LS-GI/408d | 280mm long |

Mayo Scissors



Straight


|            |            |
|------------|------------|
| LSGI/409   | 140mm long |
| LS-GI/409a | 160mm long |
| LS-GI/411  | 190mm long |
| LS-GI/411a | 230mm long |

Curved

|            |            |
|------------|------------|
| LS-GI/410  | 140mm long |
| LS-GI/410a | 160mm long |
| LS-GI/412  | 190mm long |
| LS-GI/412a | 230mm long |


Fergusson Scissors

Angled on flat



|           |            |
|-----------|------------|
| LS-GI/148 | 180mm long |
|-----------|------------|


Paynes Suture Scissors



|           |            |
|-----------|------------|
| LS-GI/152 | 100mm long |
|-----------|------------|

Ligature Stitch Scissors

Fine pattern



|           |            |
|-----------|------------|
| LS-GI/153 | 115mm long |
|-----------|------------|


Lister Bandage Scissors



|           |            |
|-----------|------------|
| LSGI/154  | 147mm long |
| LS-GI/155 | 178mm long |

McIndoe Scissors


178 mm long



|             |                  |
|-------------|------------------|
| LSGI/149    | Straight         |
| LS-GI/150   | Curved on flat   |
| LS-GI/150sc | Supercut, curved |




Heath Stitch Scissors



|           |            |
|-----------|------------|
| LS-GI/151 | 152mm long |
|-----------|------------|


Tubing Scissors



|           |            |
|-----------|------------|
| LS-GI/147 | 130mm long |
|-----------|------------|

Lister Bandage Scissors


Angled, with probe end



|           |            |
|-----------|------------|
| LS-GI/361 | 147mm long |
| LS-GI/362 | 178mm long |


Thomas Strabismus Scissors

114 mm long



|           |          |
|-----------|----------|
| LSOP/030  | Straight |
| LS-OP/029 | Curved   |

“Careless” Scissors



|           |            |
|-----------|------------|
| LS-GI/474 | 165mm long |
|-----------|------------|


Utility Scissors



|           |            |
|-----------|------------|
| LS-GI/466 | 190mm long |
|-----------|------------|

Dressing Scissors

Tungsten carbide hard insert /blunt



Straight


|           |            |
|-----------|------------|
| LS-GI/430 | 130mm long |
| LS-GI/432 | 145mm long |
| LS-GI/434 | 165mm long |
| LS-GI/436 | 185mm long |

Curved

|           |            |
|-----------|------------|
| LSGI/431  | 130mm long |
| LS-GI/433 | 145mm long |
| LS-GI/435 | 165mm long |

Dressing Scissors

Tungsten carbide hard insert, blunt/sharp



Straight


|           |            |
|-----------|------------|
| LSGI/423  | 130mm long |
| LS-GI/425 | 145mm long |
| LS-GI/427 | 165mm long |
| LS-GI/429 | 185mm long |

Curved

|           |            |
|-----------|------------|
| LS-GI/424 | 130mm long |
| LS-GI/426 | 145mm long |
| LS-GI/428 | 165mm long |

Dressing Scissors

Tungsten carbide hard inserts/sharp



Straight


|           |            |
|-----------|------------|
| LS-GI/437 | 130mm long |
| LS-GI/439 | 145mm long |
| LS-GI/441 | 165mm long |
| LS-GI/443 | 185mm long |

Curved

|           |            |
|-----------|------------|
| LS-GI/438 | 130mm long |
| LS-GI/440 | 145mm long |
| LS-GI/442 | 165mm long |

Baby Metzenbaum Scissors


Tungsten carbide, hard insert, 145mm long



|            |          |
|------------|----------|
| LSGI/445a  | Straight |
| LS-GI/446a | Curved   |

Mayo Dressing Scissors

Tungsten carbide hard insert



Straight


|            |            |
|------------|------------|
| LS-GI/419  | 140mm long |
| LS-GI/419a | 160mm long |
| LS-GI/421  | 190mm long |
| LS-GI/421a | 230mm long |

Curved

|            |            |
|------------|------------|
| LS-GI/420  | 140mm long |
| LS-GI/420a | 160mm long |
| LS-GI/422  | 190mm long |
| LS-GI/422a | 230mm long |

Mayo-Harrington Dressing Scissors

Tungsten carbide hard insert



230mm long


|           |          |
|-----------|----------|
| LSGI/417  | Straight |
| LS-GI/418 | Curved   |

200mm long

|            |          |
|------------|----------|
| LS-GI/162a | Straight |
| LS-GI/413b | Curved   |

Metzenbaum Scissors

Tungsten carbide hard insert



230mm long


|            |            |
|------------|------------|
| LS-GI/445  | 145mm long |
| LS-GI/448  | 185mm long |
| LS-GI/471a | 230mm long |
| LS-GI/471b | 250mm long |
| LS-GI/471c | 280mm long |

200mm long

|            |            |
|------------|------------|
| LS-GI/446  | 145mm long |
| LS-GI/449  | 185mm long |
| LS-GI/449a | 200mm long |
| LS-GI/449b | 230mm long |
| LS-GI/449c | 250mm long |
| LS-GI/449d | 280mm long |

Ligature Scissors


Tungsten carbide hard insert, curved



|           |            |
|-----------|------------|
| LS-GI/444 | 145mm long |
| LS-GI/447 | 185mm long |

Mayo Harrington Scissors

Supercut, tungsten carbide



230mm long


|           |          |
|-----------|----------|
| LS-GI/162 | Straight |
| LS-GI/413 | Curved   |

200mm long

|            |          |
|------------|----------|
| LS-GI/162b | Straight |
| LS-GI/413b | Curved   |

Mayo Scissors

Supercut, tungsten carbide



Straight


|            |            |
|------------|------------|
| LS-GI/163  | 140mm long |
| LS-GI/163a | 160mm long |
| LS-GI/165  | 190mm long |
| LS-GI/165a | 230mm long |

Curved

|            |            |
|------------|------------|
| LS-GI/164  | 140mm long |
| LS-GI/164a | 160mm long |
| LS-GI/166  | 190mm long |
| LS-GI/166a | 230mm long |

Metzenbaum Dissecting Scissors

Supercut, tungsten carbide



Straight


|              |            |
|--------------|------------|
| LSGI/188tc   | 185mm long |
| LS-GI/188tca | 200mm long |
| LS-GI/188tcb | 230mm long |
| LS-GI/188tcc | 250mm long |
| LS-GI/188tcd | 280mm long |

Curved

|              |         |
|--------------|---------|
| LSGI/191tc   | 1 47m m |
| LS-GI/189tc  | 185mm   |
| LS-GI/186tca | 200mm   |
| LS-GI/188tcb | 230mm   |
| LS-GI/188tcc | 250mm   |
| LS-GI/188tcd | 280mm   |

Dressing Scissors

Supercut, tungsten carbide



Blunt/Blunt Straight

|           |            |
|-----------|------------|
| LS-GI/174 | 130mm long |
| LS-GI/176 | 145mm long |
| LS-GI/178 | 165mm long |
| LS-GI/180 | 185mm long |

Curved

|            |            |
|------------|------------|
| LSGI/175   | 130mm long |
| LS-GI/177  | 145mm long |
| LS-GI/179  | 165mm long |
| LS-GI/180a | 185mm long |

Blunt/Sharp Straight

|           |            |
|-----------|------------|
| LSGI/167  | 130mm long |
| LS-GI/169 | 145mm long |
| LS-GI/171 | 165mm long |
| LS-GI/173 | 185mm long |

Curved

|            |            |
|------------|------------|
| LSGI/168   | 130mm long |
| LS-GI/170  | 145mm long |
| LS-GI/172  | 165mm long |
| LS-GI/173a | 185mm long |

Sharp/Sharp Straight

|           |            |
|-----------|------------|
| LSGI/181  | 130mm long |
| LS-GI/183 | 145mm long |
| LS-GI/185 | 165mm long |
| LS-GI/187 | 185mm long |

Curved

|            |            |
|------------|------------|
| LS-GI/182  | 130mm long |
| LS-GI/184  | 145mm long |
| LS-GI/186  | 165mm long |
| LS-GI/187a | 185mm long |

Dressing Scissors

Supercut, blunt/blunt




145mm long

|           |          |
|-----------|----------|
| LS-GI/160 | Straight |
| LS-GI/161 | Curved   |

Dressing Scissors

Supercut, blunt/sharp




145mm long

|           |          |
|-----------|----------|
| LSGI/156  | Straight |
| LS-GI/157 | Curved   |

Dressing Scissors

Supercut,sharp/sharp




145mm long

|           |          |
|-----------|----------|
| LSGI/158  | Straight |
| LS-GI/159 | Curved   |

Mayo Scissors

Supercut



Straight


|            |            |
|------------|------------|
| LSGI/166a  | 140mm long |
| LS-GI/166e | 160mm long |
| LS-GI/166b | 170mm long |

Curved

|            |            |
|------------|------------|
| LS-GI/166c | 140mm long |
| LS-GI/166f | 160mm long |
| LS-GI/166d | 170mm long |

Boyd Scissors

Supercut, blunt/blunt



|        |            |
|--------|------------|
| LS 047 | 180mm long |
|--------|------------|

Metzenbaum Dissecting Scissors

Supercut



Straight (length)

|            |       |
|------------|-------|
| LS-GI/190  | 147mm |
| LS-GI/188  | 185mm |
| LS-GI/188a | 200mm |
| LS-GI/188b | 230mm |
| LS-GI/188c | 250mm |
| LS-GI/188d | 280mm |

Curved (length)

|            |       |
|------------|-------|
| LS-GI/191  | 147mm |
| LS-GI/189  | 185mm |
| LS-GI/186a | 200mm |
| LS-GI/186b | 230mm |
| LS-GI/186c | 250mm |
| LS-GI/186d | 280mm |

Stevens Tenotomy Scissors

Supercut



Curved/blunt

|                |            |
|----------------|------------|
| LS-OP/424      | 114mm long |
| S C- O P/424 a | 125mm long |
| LS-OP/424gb    | 135mm long |

Curved/sharp

|            |            |
|------------|------------|
| LSOP/425   | 114mm long |
| LS-OP/425a | 140mm long |
| LS-OP/425b | 160mm long |

|            |                            |
|------------|----------------------------|
| LSOP/426   | Straight/blunt, 114mm long |
| LS-OP/427  | Straight/sharp, 114mm long |
| LS-OP/427c | Straight/sharp, 180mm long |

With cutout, 110mm long

|                |                |
|----------------|----------------|
| S C- O P/42 7a | Straight/blunt |
| LS-OP/427b     | Curved/blunt   |

Strabismus Scissors

Supercut




114mm long

|            |          |
|------------|----------|
| LS-GI/475  | Straight |
| LS-GI/475a | Curved   |

Iris Scissors

Supercut




114mm long

|                |          |
|----------------|----------|
| LS-OP/476      | Straight |
| S C- O P/476 a | Curved   |

Bruce Clarke Needle Holder


Screw joint



|           |            |
|-----------|------------|
| LS-GI/205 | 127mm long |
|-----------|------------|

Crile-Wood Needle Holder


Straight, baby, fine



|           |            |
|-----------|------------|
| LS-GI/199 | 150mm long |
|-----------|------------|

Kilner Needle Holder


Box joint, cross-serrated jaw



|           |            |
|-----------|------------|
| LS-GI/200 | 165mm long |
| LS-GI/201 | 133mm long |

Lawrence Needle Holder


Fine jaw, 150mm long



|            |                      |
|------------|----------------------|
| LS-PL/012b | Tungsten carbide jaw |
|------------|----------------------|

Needle Holder


box joint, lightweight,  
with curved plain jaws and one offset bow for thumb



|           |            |
|-----------|------------|
| LS-GI/209 | 152mm long |
|-----------|------------|

MacPhail Needle Holder


Screw joint, with one copper jaw



|           |            |
|-----------|------------|
| LS-GI/206 | 159mm long |
| LS-GI/207 | 178mm long |
| LS-GI/208 | 209mm long |

Mathieu Needle Holder


Screw joint



|           |            |
|-----------|------------|
| LS-GI/210 | 172mm long |
| LS-GI/211 | 209mm long |
| LS-GI/212 | 243mm long |


Mayo Needle Holder

Box joint, straight



Wide cross serrated jaws

|           |            |
|-----------|------------|
| LSGI/224  | 159mm long |
| LS-GI/398 | 203mm long |




Narrow cross serrated jaws

|           |            |
|-----------|------------|
| LS-GI/225 | 159mm long |
|-----------|------------|

Naughton Morgan Needle Holder

Box joints, cross serrated jaws




|           |            |
|-----------|------------|
| LS-GI/399 | 185mm long |
| LS-GI/400 | 254mm long |

Tungsten carbide

|           |            |
|-----------|------------|
| LS-GI/401 | 185mm long |
| LS-GI/402 | 254mm long |

Rochester-Mayo Needle Holder


Screw joint



|           |            |
|-----------|------------|
| LS-GI/213 | 178mm long |
|-----------|------------|

Berry Needle Holder/Wire Twister


Tungsten carbide jaws



|           |            |
|-----------|------------|
| LS-GI/467 | 200mm long |
|-----------|------------|

Crile-Wood Needle Holder


Tungsten carbide insert, box joint, satin finish, straight



|           |            |
|-----------|------------|
| LSGI/217  | 159mm long |
| LS-GI/218 | 178mm long |
| LS-GI/198 | 150mm long |

Derf Needle Holder


Tungsten carbide insert



|                |            |
|----------------|------------|
| LS-GI/BUPHF002 | 115mm long |
|----------------|------------|

Finochietto Needle Holder


Tungsten carbide jaws



|           |            |
|-----------|------------|
| LSGI/222  | 215mm long |
| LS-GI/223 | 265mm long |

Halsey Needle Holder


Tungsten carbide insert, box joint, satin finish, straight



|           |            |
|-----------|------------|
| LS-GI/220 | 127mm long |
|-----------|------------|

Kilner Needle Holder


Tungsten carbide insert



|           |            |
|-----------|------------|
| LS-GI/221 | 140mm long |
|-----------|------------|

Masing Needle Holder


Tungsten carbide insert, one angled bow



|            |            |
|------------|------------|
| LS-GI/009R | 130mm long |
|------------|------------|

Nievert Needle Holder


Tungsten carbide insert, box joint, satin finish, straight



|           |            |
|-----------|------------|
| LS-GI/219 | 127mm long |
|-----------|------------|

DeBakey Needle Holder


250mm long



|              |  |
|--------------|--|
| THO/LS-5642a |  |
|--------------|--|

Ryder Needle Holder

Tungsten carbide insert, box joint, satin finish, straight



0.2mm jaw pitch


|                |            |
|----------------|------------|
| LS-GI/SH-5197a | 135mm long |
| LS-GI/SH-5197b | 155mm long |
| LS-GI/SH-5197c | 175mm long |
| LS-GI/SH-5197d | 210mm long |

0.4mm jaw pitch

|                |            |
|----------------|------------|
| LS-GI/SH-5198a | 135mm long |
| LS-GI/SH-5198b | 155mm long |
| LS-GI/SH-5198c | 175mm long |
| LS-GI/SH-5198d | 210mm long |

Sarot Needle Holder


Diamond jaw, tungsten carbide, gold plated bows



|           |            |
|-----------|------------|
| LS-GI/451 | 181mm long |
|-----------|------------|

Mayo-Hegar Needle Holder

Tungstencarbideinsert,boxjoint,satinfinish, straight




|           |            |
|-----------|------------|
| LS-GI/472 | 130mm long |
| LS-GI/214 | 159mm long |
| LS-GI/215 | 178mm long |
| LS-GI/216 | 209mm long |

LONGER SIZES AVAILABLE

Olsen-Hegar Scissor/Needle Holder


With tungsten carbide insert



|           |            |
|-----------|------------|
| LSGI/229  | 147mm long |
| LS-GI/230 | 165mm long |

Gillies Scissor/Needle Holder


Right-hand, with tungsten carbide insert



|           |            |
|-----------|------------|
| LS-GI/235 | 158mm long |
|-----------|------------|

Gillies Scissor/Needle Holder


158 mm long



|           |            |
|-----------|------------|
| LS-GI/231 | Right-hand |
| LS-GI/232 | Left-hand  |

Olsen-Hegar Scissor/Needle Holder


Screw joint



|           |            |
|-----------|------------|
| LSGI/226  | 140mm long |
| LS-GI/227 | 165mm long |
| LS-GI/228 | 190mm long |

Gillies Scissor/Needle Holder


Heavy pattern, 158 mm long



|           |            |
|-----------|------------|
| LSGI/233  | Right-hand |
| LS-GI/234 | Left-hand  |

Foster Gillies Scissor/Needle Holder


Fine pattern



|           |            |
|-----------|------------|
| LS-GI/236 | 127mm long |
|-----------|------------|

Spring Clip


Serrated pattern, stainless steel, for fixing slings to leg splints, also for towel clipping



|           |           |
|-----------|-----------|
| LS-GI/344 | 63mm wide |
|-----------|-----------|

Spring Bag Clip


Serrated



|           |           |
|-----------|-----------|
| LS-GI/345 | 26mm wide |
|-----------|-----------|

Jones Towel Clip


Superior quality spring forcep type



|           |           |
|-----------|-----------|
| LS-GI/346 | 89mm long |
|-----------|-----------|

Shardles Towel Clip


Spring curved, cross action



|           |           |
|-----------|-----------|
| LS-GI/347 | 77mm long |
|-----------|-----------|

Backhaus Towel Clip


Superior quality, box joint



|            |            |
|------------|------------|
| LSGI/348   | 89mm long  |
| LS-GI/348a | 102mm long |
| LS-GI/349  | 130mm long |

Moynihan Towel Clip


Screw joint, 2-2teeth, for sides of incision



|           |            |
|-----------|------------|
| LS-GI/470 | 191mm long |
|-----------|------------|

Ball and Socket Towel Forceps

Box joint




|           |            |
|-----------|------------|
| LS-GI/350 | 89mm long  |
| LS-GI/351 | 115mm long |
| LS-GI/352 | 140mm long |

Goligher Retractor

Complete with 2 centre blades:  
1 x 90mm wide x 50mm deep  
1 x 90mm wide x 90mm deep

and 2 pairs of lateral blades:  
1 pair 25mm wide x 45mm deep  
1 pair 50mm wide x 90mm deep



|           |
|-----------|
| LS-GI/305 |
|-----------|


Balfour Abdominal Retractor

Trivalve, frame 203mm

**Child**  
102mm opening, 51mm deep lateral blades,  
38mm deep centre blade

**Adult**  
178mm opening, 64mm deep lateral blades,

**Extralarge**  
178mm opening, 102mm deep lateral blades,  
108mm deep centre blade  
44mm deep centre blade




|           |             |
|-----------|-------------|
| LSGI/373  | Child       |
| LS-GI/374 | Adult       |
| LS-GI/375 | Extra large |

Miles Retractor

Complete with 2 centre blades:  
1 x 90mm wide x 50mm deep  
1 x 90mm wide x 75mm deep

and 2 pairs of lateral blades:  
1 pair 80mm wide x 35mm deep  
1 pair 80mm wide x 75mm deep




|           |            |
|-----------|------------|
| LS-GI/376 | 340mm long |
|-----------|------------|

OTHER BLADE SIZES AVAILABLE

Balfour Abdominal Retractor


Bivale, baby



|           |                 |
|-----------|-----------------|
| LS-GI/372 | Fully open 89mm |
|-----------|-----------------|

Berkeley Abdominal Retractor


Parallel action



|           |                                     |
|-----------|-------------------------------------|
| LS-GI/461 | Blades 35mm x 60mm deep, 215mm long |
|-----------|-------------------------------------|

Turner Warwick Universal Self Retaining Retractor

|             |                                                |
|-------------|------------------------------------------------|
| LS-GI/520   | With frames A/B/C/D and blades A/B/C/D/E/F/K/L |
| LS-GI/521-A | Circular frame A only 178mm                    |
| LS-GI/521-B | Circular frame B only 229m                     |
| LS-GI/521-C | Circular frame C only 305mm                    |
| LS-GI/521-D | Oval frame D only 475mm x 330mm                |
| LS-GI/521-E | Oval frame E only 370mm x 300mm                |



Turner Warwick Needleholder with Tungsten Carbide Jaws

Curvedshankswith one short bow


|          |            |
|----------|------------|
| SH-5835  | 230mm long |
| SH-5835A | 200mm long |

Blade A




|             |                          |
|-------------|--------------------------|
| LS-GI/522-A | Blade only - 19mm x 76mm |
|-------------|--------------------------|

Blade B




|             |                           |
|-------------|---------------------------|
| LS-GI/522-B | Blade only - 38mm x 102mm |
|-------------|---------------------------|

Blade C




|             |                          |
|-------------|--------------------------|
| LS-GI/522-C | Blade only - 57mm x 89mm |
|-------------|--------------------------|

Blade D




|             |                           |
|-------------|---------------------------|
| LS-GI/522-D | Blade only - 57mm x 114mm |
|-------------|---------------------------|

Blade E




|             |                                                       |
|-------------|-------------------------------------------------------|
| LS-GI/522-E | Blade only - curved stainless steel blade 51mm x 38mm |
|-------------|-------------------------------------------------------|

Blade F




|             |                                                       |
|-------------|-------------------------------------------------------|
| LS-GI/522-F | Blade only - curved stainless steel blade 76mm x 60mm |
|-------------|-------------------------------------------------------|

Blade G




|             |                                                     |
|-------------|-----------------------------------------------------|
| LS-GI/522-G | Blade only - flat stainless steel blade 89mm x 60mm |
|-------------|-----------------------------------------------------|

Blade J




|             |                                     |
|-------------|-------------------------------------|
| LS-GI/522-J | Blade only - single hook, 13mm deep |
|-------------|-------------------------------------|

Blade K



|             |                                         |
|-------------|-----------------------------------------|
| LS-GI/522-K | Blade only - flat with clip 55mm x 90mm |
|-------------|-----------------------------------------|


Blade L



|             |                                          |
|-------------|------------------------------------------|
| LS-GI/522-L | Blade only - flat with clip 55mm x 120mm |
|-------------|------------------------------------------|

Collin Retractor


110mm fully open, 235mm long



|            |                                                                     |
|------------|---------------------------------------------------------------------|
| LS-GI/322  | Complete with centre blade 55x40mm and one pair of blades 38mmx60mm |
| LS-GI/322a | Without centre blade                                                |


Gelpi Retractor

Selfretaining



|            |           |
|------------|-----------|
| LSGY/013a  | 140mmlong |
| LS-GY/013b | 178mmlong |


Pozzi Abdominal Retractors



|           |                                             |
|-----------|---------------------------------------------|
| LS-GI/370 | Swivel blades: 38x 60mm, bivalve,254mm long |
|-----------|---------------------------------------------|


Joll Thyroid Retractor

Self retaining, with lock-innate



|           |           |
|-----------|-----------|
| LS-GI/363 | 152mmlong |
|-----------|-----------|


Gossett Abdominal Retractor



|           |                           |
|-----------|---------------------------|
| LS-GI/303 | Large,fully open185mmlong |
| LS-GI/304 | Small,fully open95mmlong  |

Bowlby Dilator


Curved



|           |           |
|-----------|-----------|
| LS-GI/495 | 145mmlong |
|-----------|-----------|

Irwin Retractor


Self retaining, with cam rack, 4 x 4 teeth



|           |           |
|-----------|-----------|
| LS-GI/324 | 178mmlong |
|-----------|-----------|

Norfolk and Norwich Retractor


4 x 5 teeth 50mm blades x 223mm long



|           |               |
|-----------|---------------|
| LS-GI/325 | Standard rack |
|-----------|---------------|

Travers Retractor


Self retaining with internal rack, 210mm long



|           |                        |
|-----------|------------------------|
| LS-GI/323 | 4x5teeth,95mmfullyopen |
|-----------|------------------------|

West Retractor


Self retaining, 140mm long



|             |                                     |
|-------------|-------------------------------------|
| LS-GI/326   | 3x4teeth,60mmfullyopen              |
| LS-GI/326cr | 3x4teeth,60mmfullyopen,withcam rack |

Weitlander Retractor


3-4 teeth



| 135mm long |                |
|------------|----------------|
| LSGI/496   | Blunt          |
| LS-GI/496a | Sharp          |
| LS-GI/496b | Sharp,2-3teeth |
| LS-GI/496c | Blunt,2-3teeth |

| 165mm long |       |
|------------|-------|
| LS-GI/497  | Blunt |
| LS-GI/497a | Sharp |

Kocher Retractor



230mm long


|           |                 |
|-----------|-----------------|
| LS-GI/378 | 36 x 20mm blade |
|-----------|-----------------|

250mm long

|            |                 |
|------------|-----------------|
| LS-GI/378a | 67 x 25mm blade |
| LS-GI/378b | 80 x 42mm blade |
| LS-GI/378c | 80 x 50mm blade |
| LS-GI/378d | 80 x 60mm blade |

Canny Ryall Navy Retractor


Double-ended, with one single hook blade 5mm wide, and one fenestrated blade



|            |                                         |
|------------|-----------------------------------------|
| LSGI/343   | Fenestrated blade 16mm wide, 190mm long |
| LS-GI/343a | Fenestrated blade 19mm wide, 190mm long |
| LS-GI/343b | Fenestrated blade 25mm wide, 220mm long |
| LS-GI/343c | Fenestrated blade 30mm wide, 220mm long |

Morris Retractor


Single blade, 229mm long



|           |           |
|-----------|-----------|
| LS-GI/330 | 38 x 51mm |
| LS-GI/331 | 51 x 51mm |
| LS-GI/332 | 64 x 51mm |

Morris Retractor


Double ended, 216 mm long



|           |                         |
|-----------|-------------------------|
| LSGI/327  | 32 x 51mm and 38 x 51mm |
| LS-GI/328 | 44 x 51mm and 51 x 51mm |
| LS-GI/329 | 64 x 51mm and 38 x 51mm |

Kocher Langenbeck Retractor

215 mm long




|             |                 |
|-------------|-----------------|
| LSGI/336sp  | 25 x 6mm blade  |
| LS-GI/379a  | 35 x 8mm blade  |
| LS-GI/379b  | 35 x 11mm blade |
| LS-GI/379c  | 35 x 15mm blade |
| LS-GI/357sp | 41 x 11mm blade |


|                |                 |
|----------------|-----------------|
| LSTHO/cutsp015 | 55 x 11mm blade |
| LS-GI/379d     | 70 x 15mm blade |
| LS-GI/379e     | 80 x 12mm blade |
| LS-GI/379f     | 80 x 16mm blade |

Canny Ryall Army Retractor


Double prong, 13mm wide, 165mm long



|            |                  |
|------------|------------------|
| LS-GI/342  | Large, 19mm wide |
| LS-GI/342a | Small, 11mm wide |



Moyniham Retractor



|           |                                           |
|-----------|-------------------------------------------|
| LS-GI/462 | Single ended, blade 13 x 22mm, 215mm long |
|-----------|-------------------------------------------|


Deaver Abdominal Retractors



|           |                                |
|-----------|--------------------------------|
| LS-GI/306 | Set of 5, 315mm overall length |
| LS-GI/307 | 25mm wide                      |
| LS-GI/308 | 25mm wide                      |
| LS-GI/309 | 38mm wide                      |
| LS-GI/310 | 51mm wide                      |
| LS-GI/311 | 76mm wide                      |

Langenbeck Retractor


Single ended



|             |                   |
|-------------|-------------------|
| LSGI/336    | Blade 22mm x 6mm  |
| LS-GI/337   | Blade 45mm x 13mm |
| LS-GI/337a  | Blade 65mm x 13mm |
| LS-GI/338a  | Blade 75mm x 25mm |
| LS-GI/338a1 | Blade 65mm x 25mm |


Czerny Retractor

Double ended



|           |            |
|-----------|------------|
| LS-GI/341 | 171mm long |
|-----------|------------|

Elliot (Modified Deaver) Retractor



|            |                  |
|------------|------------------|
| LS-GI/308M | With flat handle |
|------------|------------------|

Dyball and Lock Retractor


For liver




|           |                        |
|-----------|------------------------|
| LS-GI/335 | 64 x 144mm, 241mm long |
|-----------|------------------------|

OTHER SIZES AVAILABLE - PRICE ON APPLICATION

Landon Retractor




|           |       |
|-----------|-------|
| LS-GI/494 | 230mm |
|-----------|-------|



Lane Retractor


Double ended blade, 229mm long



|           |                  |
|-----------|------------------|
| LS-GI/339 | 13mm x 19mm wide |
| LS-GI/340 | 25mm x 32mm wide |

Deaver Abdominal Retractors

With round handles



|           |                                |
|-----------|--------------------------------|
| LS-GI/312 | Set of 5, 305mm overall length |
| LS-GI/313 | 25mm wide                      |
| LS-GI/314 | 25mm wide                      |
| LS-GI/315 | 38mm wide                      |
| LS-GI/316 | 51mm wide                      |
| LS-GI/317 | 76mm wide                      |

Kelly Retractor Blade


267 mm long



|           |            |
|-----------|------------|
| LSGI/333  | 38 x 191mm |
| LS-GI/334 | 57 x 191mm |

Durham Retractor

With fenestrated blade, 216mm long



|           |           |
|-----------|-----------|
| LS-GI/300 | 10mm wide |
| LS-GI/301 | 14mm wide |
| LS-GI/302 | 20mm wide |

Ballance (Volkman) Retractor

203mm long




Sharp

|            |         |
|------------|---------|
| LSGI/463   | 1 prong |
| LS-GI/463a | 2 prong |
| LS-GI/463b | 3 prong |
| LS-GI/463c | 4 prong |
| LS-GI/463d | 6 prong |

Blunt

|            |         |
|------------|---------|
| LSGI/463e  | 1 prong |
| LS-GI/463f | 2 prong |
| LS-GI/463g | 3 prong |
| LS-GI/463h | 4 prong |
| LS-GI/463i | 6 prong |


Volkman Retractor



|            |                   |
|------------|-------------------|
| LS-OR/400  | Sharp, 220mm long |
| LS-OR/400a | Blunt, 220mm long |
| LS-OR/400b | Sharp, 245mm long |


Kilner Retractor

Double ended, cat's paw



|           |            |
|-----------|------------|
| LS-PL/056 | 152mm long |
|-----------|------------|


Kocher Retractor



|             |                   |
|-------------|-------------------|
| LS-OR/N036c | Sharp, 205mm long |
| LS-OR/N036  | Sharp, 220mm long |
| LS-OR/N036b | Blunt, 220mm long |

Senn-Miller Retractor

170mm long



|            |
|------------|
| LS-SH/5881 |
|------------|

Single Hook Retractor

Sharp



|              |                 |
|--------------|-----------------|
| LS-GI/SH5671 | 152mm long      |
| LS-PL/065b   | Half size, fine |

Gillies Skin Hook

Solid handle, 190mm long




|            |                 |
|------------|-----------------|
| LS-PL/065a | Standard, fine  |
| LS-PL/065b | Half size, fine |

MacDonald Dissector Double Ended Probe, Sinus Silver



|           |            |
|-----------|------------|
| LS-GI/285 | 191mm long |
|-----------|------------|


Brodie Director Probe Ended, Silver Plated



|           |            |
|-----------|------------|
| LS-GI/286 | 165mm long |
|-----------|------------|

Watson Cheyne Dissector

Fine, double-ended, with probe end



|            |            |
|------------|------------|
| LS-GI/287  | 127mm long |
| LS-GI/288  | 191mm long |
| LS-GI/288a | 230mm long |
| LS-GI/288b | 255mm long |

Probe, Malleable Silver Plated




|            |            |
|------------|------------|
| LS-GI/289a | 127mm long |
| LS-GI/289  | 152mm long |
| LS-GI/290  | 178mm long |
| LS-GI/291  | 203mm long |
| LS-GI/292  | 255mm long |
| LS-GI/292a | 305mm long |



|           |            |
|-----------|------------|
| LS-GI/464 | 152mm long |
|-----------|------------|

Probe, Fish Tail Silver Plated



|           |            |
|-----------|------------|
| LSGI/293  | 127mm long |
| LS-GI/294 | 152mm long |
| LS-GI/295 | 203mm long |
| LS-GI/296 | 254mm long |

Jobson Horne Probe Dressing (Aural)

With one flat tapered end, serrated edge and one small ring end, stainless steel



|           |            |
|-----------|------------|
| LS-GI/297 | 140mm long |
| LS-GI/298 | 178mm long |


Aneurysm Needle

Small, for infusions




|            |            |
|------------|------------|
| LS-GI/472a | 146mm long |
|------------|------------|

Syme Aneurysm Needle




|           |            |
|-----------|------------|
| LS-GI/473 | 165mm long |
|-----------|------------|

Ring Cutter



|            |                                |
|------------|--------------------------------|
| LS-GI/452  | Circular saw blade, 165mm long |
| LS-GI/452a | Spare blade                    |


Mayo Instrument Pin



|           |            |
|-----------|------------|
| LS-GI/355 | 102mm long |
| LS-GI/356 | 152mm long |

Baron Scalpel Handle

130mm long



|           |      |
|-----------|------|
| LS-PM/003 | No 3 |
| LS-PM/004 | No 4 |

Scalpel Handle


127mm long



|            |        |
|------------|--------|
| LS-PM/007  | No 3   |
| LS-PM/008  | No 4   |
| LS-PM/008a | No 4 L |
| LS-PM/001a | No 7   |

Kidney Dish

Stainless steel



|           |       |
|-----------|-------|
| LS-GI/380 | 152mm |
| LS-GI/381 | 203mm |
| LS-GI/382 | 254mm |
| LS-GI/383 | 305mm |

Lotion Bowl


Without lid, stainless steel



|           |       |
|-----------|-------|
| LS-GI/384 | 102mm |
| LS-GI/385 | 152mm |
| LS-GI/386 | 203mm |
| LS-GI/387 | 254mm |
| LS-GI/388 | 305mm |

Instrument Tray

Stainless steel



|           |                               |
|-----------|-------------------------------|
| LS-GI/389 | 254 x 203 x 50mm              |
| LS-GI/390 | Lid only for 254 x 203 x 50mm |
| LS-GI/391 | 305 x 254 x 50mm              |
| LS-GI/392 | Lid only for 305 x 254 x 50mm |
| LS-GI/393 | 356 x 305 x 50mm              |
| LS-GI/394 | Lid only for 356 x 305 x 50mm |

Cheate Jar



|           |                               |
|-----------|-------------------------------|
| LS-GI/397 | Stainless steel, 203mm x 50mm |
|-----------|-------------------------------|

Measuring Jug




|           |                  |
|-----------|------------------|
| LS-GI/395 | 1 litre capacity |
|-----------|------------------|

Petri Dish



|              |                                  |
|--------------|----------------------------------|
| LS-SH/002149 | 100mm, stainless steel, with lid |
|--------------|----------------------------------|

Gallipot



|           |                 |
|-----------|-----------------|
| LS-GI/396 | Stainless s63mm |
|-----------|-----------------|

Sterilising Baskets

DIN Basket

Wire mesh base with perforated sides and two carrying handles

|             |                  |             |                  |
|-------------|------------------|-------------|------------------|
| LS-GI/CO006 | 230 x 130 x 60mm | LS-GI/CO013 | 480 x 250 x 70mm |
| LS-GI/CO007 | 240 x 130 x 50mm | LS-GI/CO016 | 480 x 400 x 60mm |
| LS-GI/CO008 | 240 x 130 x 60mm | LS-GI/CO017 | 480 x 480 x 50mm |
| LS-GI/CO009 | 240 x 250 x 50mm | LS-GI/CO018 | 480 x 480 x 60mm |
| LS-GI/CO010 | 240 x 250 x 60mm | LS-GI/CO019 | 500 x 300 x 50mm |
| LS-GI/CO014 | 300 x 300 x 50mm | LS-GI/CO020 | 550 x 430 x 50mm |
| LS-GI/CO015 | 350 x 350 x 60mm | LS-GI/CO021 | 600 x 300 x 60mm |
| LS-GI/CO011 | 480 x 250 x 50mm |             |                  |
| LS-GI/CO012 | 480 x 250 x 60mm |             |                  |

DIN Basket

Wire mesh basket with two carrying handles

|             |                  |             |                   |
|-------------|------------------|-------------|-------------------|
| LS-GI/CO022 | 240 x 250 x 50mm | LS-GI/CO026 | 480 x 250 x 60mm  |
| LS-GI/CO023 | 240 x 250 x 60mm | LS-GI/CO027 | 480 x 250 x 70mm  |
| LS-GI/CO024 | 240 x 250 x 70mm | LS-GI/CO028 | 480 x 250 x 100mm |
| LS-GI/CO031 | 300 x 300 x 60mm | LS-GI/CO029 | 480 x 480 x 50mm  |
| LS-GI/CO025 | 480 x 250 x 50mm | LS-GI/CO030 | 480 x 480 x 70mm  |


DIN Basket Lids

For perforated and wire mesh baskets

|             |             |
|-------------|-------------|
| LS-GI/CO032 | 240 x 250mm |
| LS-GI/CO033 | 300 x 300mm |
| LS-GI/CO034 | 480 x 250mm |
| LS-GI/CO035 | 480 x 480mm |

ENT Basket

Wire mesh base with perforated sides and two carrying handles, supplied with lid



|             |                   |
|-------------|-------------------|
| LS-GI/CO036 | 480 x 250 x 160mm |
| LS-GI/CO037 | Fixing clamps     |

Arthroscope Basket

Wire mesh base with two fixing clamps



|             |                 |
|-------------|-----------------|
| LS-GI/CO001 | 290 x 80 x 52mm |
|-------------|-----------------|

Endoscope Basket


With small lens basket

|             |                   |
|-------------|-------------------|
| LS-GI/CO002 | 6400 x 150 x 77mm |
|-------------|-------------------|

Wire mesh base with two fixing clamps

|             |                     |
|-------------|---------------------|
| LS-GI/CO003 | S - 460 x 80 x 52mm |
| LS-GI/CO004 | D - 480 x 80 x 52mm |
| LS-GI/CO005 | D - 680 x 80 x 52mm |

Sieve



|            |                |
|------------|----------------|
| LS-GI/5071 | 203mm x 50mm ø |
|------------|----------------|

Fine Mesh Basket



With hinged take off / lockable lid With separate lid

|             |                  |             |                  |
|-------------|------------------|-------------|------------------|
| LS-GI/CO041 | 120 x 155 x 40mm | LS-GI/CO038 | 120 x 155 x 40mm |
| LS-GI/CO044 | 120 x 155 x 50mm | LS-GI/CO039 | 235 x 155 x 40mm |
| LS-GI/CO042 | 235 x 155 x 40mm | LS-GI/CO040 | 375 x 260 x 40mm |
| LS-GI/CO045 | 235 x 155 x 50mm |             |                  |
| LS-GI/CO043 | 375 x 260 x 40mm |             |                  |

Contact us for full spec and to order  
info@surgicalholdings.co.uk