



PLASTIC / ORAL SURGERY



A versatile range for reconstructive, aesthetic, and maxillofacial procedures. Includes fine scissors, tissue forceps, needle holders, elevators, periosteal raspatories, and self-retaining retractors. Designed for precise soft tissue handling, delicate suturing, and controlled bone contouring. Combines balanced ergonomics, fine tips, and durable stainless steel construction to support meticulous work in both superficial and deep facial planes.

Foster Scissor/Needle Holder

Fenestrated jaw, 125mm long



LS-PL/010a	Left
------------	------

Derf Needle Holder


Tungsten carbideinsert, with two gold bows



LS-PL/013	115mm long
LS-PL/013L	125mm long

Masing Needle Holder


Tungstencarbide insert, one angled bow



LS-GI/009R	130mm long
------------	------------

Castroviejo Needle Holder

LS-OP/021	140mm long, with catch
LS-OP/020	140mm long, without catch



Diamond dust jaws

LS-THO/060	180mm long, without catch
LS-THO/061	200mm long, without catch
LS-THO/062	215mm long, without catch

ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO


Tungsten carbide jaws

LS-THO/058a	180mm long
LS-THO/058	200mm long
LS-THO/059	215mm long

ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO

Lawrence Needle Holder


Fine jaw, 150mm long



LS-PL/012b	Tungsten carbide, with two gold bows
------------	--------------------------------------

Reverdin Suture Needle


With automatic eye, 203mm long



LS-PL/013a	Fully curved
LS-PL/013b	Curved
LS-PL/013c	Slightly curved
LS-PL/013d	Straight

Mustarde Suture Holding Forceps


Convex and concave jaw



LS-PL/011	133mm long
-----------	------------


Fickling Artery Forceps

190mm long



LS-PL/044a	Serrated
LS-PL/044b	1-2 teeth


Asch Septum Forceps



LS-PL/031	235mm long
-----------	------------

Jansen-Middleton Septum Forceps


Compound action, 178mm long



LS-PL/032	5mm bite
-----------	----------

Walsham Septum Straightening Forceps


220mm long



LS-PL/030a	Straight
LS-PL/030b	Right
LS-PL/030c	Left

Bone Cutting Forceps


Straight, compound action



LS-PL/033	146mm long
-----------	------------

Ward Bone Cutting Forceps

Slightly angled onflat



LS-PL/038	178mm long
-----------	------------


Cleft Palette Raspatory



LS-PL/091	Curved downwards
LS-PL/092	Curved upwards

Cleft Palette Raspatory


Curved upwards




LS-PL/091	Left
LS-PL/092	Right

Aufricht Raspatory

9mm wide, 220mm long




LS-PL/095	Up stroke
LS-PL/095a	Down stroke



Gallagher Raspatory

190mm long



LS-PL/098

McIndoe Raspatory

Spoonend, 170mm long



LS-/PL088	Small
LS-/PL089	Medium
LS-/PL090	Large

Farabeuf Elevator

Light, 160mm long



LS-PL/024	12mm curved
LS-PL/024a	8mm curved

Durogrip Raspatory

7mm wide, 205mm long






12/14 teeth

LS-PL/096	12/14 teeth
LS-PL/096a	16/20 teeth
LS-PL/096b	24/30 teeth
LS-PL/096c	36/44 teeth
LS-PL/096d	52/62 teeth


Joseph Raspatory

Chequered, 165mm long



LS-PL/099	Fine, 7.5mm
LS-PL/099a	Coarse, 9mm

Howarth Elevator



LS-PL/025	220mm long
-----------	------------

Syme Periosteal Elevator


Double ended



LS-PL/046	170mm long
-----------	------------

Cryer Elevator


150mm long



LS-PL/069	Right
LS-PL/069a	Left

Coupland Elevator


Curved, 150mm long



LS-PL/068	3mm wide
LS-PL/068a	3.5mm wide
LS-PL/068b	4mm wide

Fomon Raspatory



6.5 and 7.5mm wide, 215mm long



LS-PL/097	Chequered
LS-PL/097a	Knotched

Lewis Raspatory

8mm wide, 175mm long




LS-PL/100	Fine teeth
LS-PL/100a	Coarse teeth

Maltz Raspatory

180mm long

LS-PL/101	Straight, 9mm
LS-PL/101b	Straight, 5mm
LS-PL/101c	Curved, 5mm




Miller Raspatory

180mm long



LS-PL/102	4.3 & 6mm
LS-PL/102a	5 & 6mm
LS-PL/102b	5 & 8mm



Gillies Dissecting Forceps


LS-PL/052	1 x 2 teeth, 152mm long
-----------	-------------------------



McIndoe Dissecting Forceps

152mm long

LS-PL/053a	1-2 teeth
LS-PL/053b	Serrated



Waugh Dissecting Forceps

Fine, serrated


LS-GI/261	152mm long
LS-GI/262	178mm long
LS-GI/263	203mm long



College Pattern Stopping and Dressing Forceps

152mm long

LS-PL/048	Serrated
-----------	----------



Mustarde Dissecting and Suture Holding Forceps

Combined double hooked jaws


LS-PL/050	156mm long
-----------	------------



Childe Approximating Forceps

With magazine for Michel clips

LS-PL/073	180mm long
-----------	------------



Mitchell Trimmer

165mm long

LS-ENT/013



Dental Pick

LS-PL/067	160mm long
-----------	------------



Gillies Skin Hook

190mm long

LS-PL/065a	Standard
LS-PL/065b	Half Size



Joseph Double Skin Hook

165mm long

LS-PL/056a




Kleinert Kutz Skin Hooks

LS-PL/109	16cm large
LS-PL/108	16cm medium
LS-PL/107	16cm small



Alms Skin Retractor


Self retaining



LS-PL/060	70mm long
LS-PL/060a	100mm long

Kilner Retractor


Double ended, cat's paw



LS-PL/056	152mm long
-----------	------------

Kilner Alae Retractor

85mm long



LS-PL/063a	10mm
LS-PL/063b	13mm

Austin Kay Retractor


85mm long



S C- P L /074	16mm wide
---------------	-----------

Bowdler Henry Retractor


170mm long



S C- P L /076	12mm wide
---------------	-----------

Tongue Retractor


190mm long



LS-PL/071	12mm wide
-----------	-----------

Cheek Retractor

140mm long



LS-PL/078	40mm wide
-----------	-----------

Ragnell Davis Retractor


140mm long



LS-PL/074a

Flap Retractor

180mm long



LS-PL/079	32mm wide
-----------	-----------

Kilner Nasal Retractor


152mm long



LS-PL/061a	7mm wide
------------	----------

Hajek Retractor


Upper lip, 165mm long



LS-PL/059	20mm wide
-----------	-----------

Freece Oral Retractor

150mm long



LS-PL/080	33mm wide
-----------	-----------

Kilner Cheek Retractor


Double ended,167mm long



LS-PL/058a	Small, 19mm and 25mm wide
LS-PL/058b	Large, 29mm and 35mm wide

Laster Cheek Retractor

110mm long



LS-PL/081	35mm wide
-----------	-----------

Ward (East Grinstead) Osteotome


Square handle, 178mm long



LS-PL/040aeg	3mm wide
LS-PL/040beg	5mm wide

Gillies (French Pattern) Osteotome

140mm long



LS-PL/026	3mm wide
LS-PL/026a	5mm wide
LS-PL/026b	7mm wide
LS-PL/026c	9mm wide
LS-PL/026g	11mm wide
LS-PL/026d	16mm wide

Jenkins Chisel


140mm long



LS-PL/027a	2mm wide
LS-PL/027b	4mm wide
LS-PL/027c	6mm wide
LS-PL/027d	8mm wide
LS-PL/027e	10mm wide

McIndoe Nasal Chisel

152mm long



LS-PL/028a	5mm wide
LS-PL/028b	7mm wide
LS-PL/028c	9mm wide
LS-PL/028d	11mm wide


Silver Nasal Chisel

182mm long



LS-PL/029a	Straight
LS-PL/029b	Curved to right
LS-PL/029c	Curved to left


Pterigoideus (Maxilla) Chisel



LS-PL/087	255mm long
-----------	------------

Tilley Septum Gouge


V shaped, straight



LS-PL/043	127mm long
-----------	------------


Fickling Dental Mallet

191mm long




LS-PL/044	25mm ø, tapering head
-----------	-----------------------

Silver Skin Graft Knife



LS-PL/015a	190mm long
------------	------------

Eckhoff Mapping Pen



LS-PL/019	145mm long
LS-PL/019a	Nibs packs of 10

Sommerland Surgical Marking Pen


Adjustable four-pointed tip



LS-PL/020	152mm long
-----------	------------


Gillies Suction Tube

Control valve, 180mm long




LS-PL/075a	2.5mm
LS-PL/075b	3mm
LS-PL/075c	3.5mm

Watson Skin Graft Knife




LS-PL/017	Right-handed
LS-PL/017L	Left-handed

Doyen Mouth Gag




LS-ENT/88d	Adult
LS-ENT/88e	Child

Fergusson Ackland Mouth Gag



LS-ENT/209	Adult, 210mm
LS-ENT/209a	Child, 155mm

Dental Cartridge Syringe




LS-PL/072	115mm long
-----------	------------

Freeman Areola Nipple Sizer


1.5mm

LS-PL/NIP25	25mm ø
LS-PL/NIP40	40mm ø
LS-PL/NIP45	45mm ø
LS-PL/NIP50	50mm ø
LS-PL/NIP55	55mm ø



Kelsey Fry Bone Awl


148mm long



LS-PL/039a	Straight
LS-PL/039b	Curved

Joseph Nasal Saw

Blade 35mm long



Horizontal shaft

LS-PL/034a	Cranked right
LS-PL/034b	Cranked left

Vertical shaft

LS-PL/034c	Cranked right
LS-PL/034d	Cranked left

Unna Comedo Expressor Double Ended Curette

2mm / 1mm holes



LS-PL/035	135mm long
-----------	------------

Lang Steel Scoop

Double ended



LS-PL/064	7mm x 4mm, 140mm long
-----------	-----------------------

McHardy Curette

136mm long



LS-PL/036	1mm, 1.5mm, and 2mm ø
-----------	-----------------------

MacEwen Curette


Combined with seeker, 165mm long



LS-PL/045a	Oval
------------	------

Kilner Scissors


Rounded blades, 115mm long



LS-PL/001a	Straight, sharp points
LS-PL/001b	Curved, blunt points

Aufricht Dissecting Scissors


Heavy pattern, triangular blades, 140mm long



LS-PL/002a	Straight
LS-PL/002d	Curved

Aufricht Dissecting Scissors


Light pattern, 150mm long



LS-PL/002b	Straight
LS-PL/002c	Curved

Scissors


Curved, sharp points



LS-PL/004	130mm long
-----------	------------

Scissors

Angular, for cutting cartilage, wire cutting notch



LS-PL/006	127mm long
-----------	------------

McIndoe Scissors


Curved on flat, round tapering ends



LS-PL/007	190mm long
-----------	------------

Stevens Tenotomy Scissors


Super cut, 114mm long



LS-OP/424	Curved/blunt
LS-OP/425	Curved/sharp
LS-OP/426	Straight/blunt
LS-OP/427	Straight/sharp

Cartilage Cutting Scissors


Finepoints, serrated edges



LS-PL/009	190mm long
-----------	------------

Gorney Supercut Scissors


Rhytidectomy scissor, tungsten carbide, one micro serrated and one razor sharp edge



LS-PL/085a	125mm long
LS-PL/085	200mm long, straight
LS-PL/085b	200mm long, curved

Fickling (Cathcart) Scissors

Angled



LS-PL/049	180mm long
-----------	------------