

PROCTOLOGY

A focused selection of instruments for diagnostic and surgical procedures of the anorectal region. Includes proctoscopes, anosopes, sigmoidoscopes, biopsy forceps, hemorrhoidal clamps, ligators, and retractors. Designed to provide optimal visualization and access within confined anatomy while ensuring atraumatic tissue handling. Crafted from premium stainless steel with smooth, polished surfaces for patient comfort and reliable, long-term performance.



Officer Biopsy Forceps

With 8mm cup jaw, to use through a sigmoidoscope



Yeoman Basket

S-Jawed to use through a sigmoidoscope




Lloyd-Davis Sigmoidoscope

Biopsy forceps



Walton Rectal Biopsy Forceps

With laterally angled shanks




Officer Dissecting Forceps




Lloyd-Davis Rectal Occlusion Forceps

With longitudinal serrations




Finch Excision Clamp

With hollowed jaws



Hays Colon Resection Clamp


Angled at 70°



Abel Abdominal-Perineal Scissors

Straight

LS-PR/016	320mm long
-----------	------------




Goligher Scissors

Straight

LS-PR/017a	300mm long
LS-PR/017b	330mm long

Curved


LS-PR/017c	300mm long
LS-PR/017d	330mm long



Lloyd-Davis Rectal Scissors


Flattened blades, blunt points
265mm long

LS-PR/018a	Straight
LS-PR/018b	Curved on flat



St Mark's Rectal Dilator

LS-PR/005a	Small, 14mm x 19mm
LS-PR/005b	Large, 16mm x 25mm



Lockhart Mummery Fistula Director

165mm long


Straight

LS-PR/006a1	With 1mm hole
-------------	---------------




90° curve

LS-PR/006c1	With 1mm hole
-------------	---------------




45° curve

LS-PR/006b	With 1mm hole
------------	---------------




Retrograde

LS-PR/006d1	With 1mm hole
-------------	---------------




Lloyd-Davis Retractor

LS-PR/013a	Short, without tapered end, 125mm x 55mm
LS-PR/013b	Long, with tapered end, 170mm x 65mm x 50mm



Park Combined Anal Retractor


Self retaining, with 2 sets of interchangeable blades and one adjustable centre blade



LS-PR/015	140mm long
-----------	------------

Park Anal Retractor


Self retaining



LS-PR/014	130mm, 95mm
-----------	-------------


Kelly Proctoscope

With obturator, lumen



LS-PR/012a	140mm long
LS-PR/012b	200mm long
LS-PR/012c	Tapered, 22mm ø, x 100mm

Eisenhammer Speculum



LS-SH/4120	Blade 70mm long x 20mm wide
LS-SH/41201	Blade 95mm long x 20mm wide